

CTS-5531

Solvent Reservoir Filter, 10µm
Comparable to: Waters®
WAT025531 and 700003615



Technical Data Sheet

1.0 Description

This 10µm solvent reservoir filter is designed to be compatible with Waters® 510, 515, 600, 610, 1515, 1525, 1525 Micro, 1525EF, 2535, 2545, 2555, 2690, 2690D, 2695, 2695D, 2700 and Alliance® HPLC systems. It is also compatible with Acquity® and Acquity UPLC® BSM UHPLC systems. It is comparable to Waters® part number WAT025531 and 700003615.

2.0 Materials Specifications

316 SS

3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.



4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at sales@sciencix.com to report the issue. We will assist you with a resolution in a timely fashion.

Please note: A short piece of 1/8" ID tubing, such as CTS-6809 (Mobile Phase Transfer Line, TFE, 3/16" OD x 1/8" ID x 25'), is needed to attach the CTS-5531 solvent filter to standard 1/8" OD solvent tubing.