

# CTS-21650

Performance Maintenance Kit, e2695  
Comparable to: Waters® 201000313



## Technical Data Sheet

### 1.0 Description

This PM kit contains all the parts necessary to keep Waters® e2695 systems working at peak performance. It is comparable to Waters® part number 201000313.

### 2.0 Specifications

#### 2.1 Components

Qty	Part Description
2	Sapphire Plunger (CTS-10476)
1	Seal Wash Plunger Seal Kit, 2/pk (CTS-10560)
1	Head Plunger Seal Kit, 2/pk (CTS-10566)
1	Wash Tube Seal Kit, 4/pk (CTS-10597)
1	Sparge Diffuser (CTS-10984)
1	Face Seal, 4/pk (CTS-10989)
5	Solvent Reservoir Filter, 10µm (CTS-5531)
1	250µL Syringe (CTS-3109)
1	2µm Precolumn Filter Insert (CTS-10919)
2	Seal Wash Face Seal (CTS-10532)
1	Check Valve Cartridge, 2/pk (CTS-10600)
1	Performance Autosampler Needle (CTS-7550)
4	Washer, Small, PTFE (CTS-7321)
1	Compression Screw, 1/16", SS (CTS-5313)
1	Needle Wash Frit (CTS-2324)
1	Gold Injector Seal, 2/pk (CTS-A047)
1	Ferrule, .061 (CTS-21654)
1	Washer, PTFE (CTS-10741-01)
2	Zip Tie-White (CTS-21432)
1	Label, Field-Service (CTS-21447)



#### 2.2 Materials

18-8 SS	HDPE
302 SS	MP35N SS
316 SS	PCTFE
Borosilicate Glass	PEEK™
Delrin®	PTFE
Elgiloy Spring	Ruby
UHMWPE	Sapphire

#### 3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

#### 4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at [sales@sciencix.com](mailto:sales@sciencix.com) to report the issue. We will assist you with a resolution in a timely fashion.