

# CTS-21693

Tungsten Lamp, i-Series

Comparable to: Shimadzu®: 228-34410-41, Old # 062-65012-04,

Old #228-34410-91



## Technical Data Sheet

---

### 1.0 Description

This lamp is designed to be compatible with Shimadzu® LC-2030 3D, LC-2040 3D detectors. It is comparable to Shimadzu® part numbers 228-34410-41, Old # 062-65012-04, Old # 228-34410-91.

### 2.0 Material Specifications

Tungsten Halogen

### 3.0 Replacement Frequency

It is recommended that the lamp be replaced every two years, or if it will not ignite, the noise or drift exceeds specifications, or if the lamp energy is below 50% of original intensity, as recorded during the initial installation.

### 4.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.



### 5.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at [sales@sciencix.com](mailto:sales@sciencix.com) to report the issue. We will assist you with a resolution in a timely fashion.