

# CTS-11093-02

Plunger, 2/pk  
Comparable to: Waters®  
700002600



## Technical Data Sheet

---

### 1.0 Description

These plungers are designed to be compatible with Waters Acquity®, Acquity H-Class Bio QSM, Acquity H-Class QSM, Acquity I-Class BSM, Acquity Isocratic Solvent Mgr, Acquity M-Class μASM, Acquity M-Class μBSM, Acquity UPLC® BSM, Acquity UPLC I2V BSM, nanoAcquity® and nanoAcquity UPLC BSM/ASM UHPLC systems. They are comparable to Waters part number 700002600.



### 2.0 Materials Specifications

Sapphire  
316 SS

### 3.0 Replacement Frequency

Replacement may be indicated by a leak at the pump head. The sapphire plunger should be inspected with a magnifier whenever the pump head is disassembled to replace the plunger seals. If any scratches are observed or a worn or flat area seen the plunger should be replaced.

### 4.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

### 5.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at [sales@sciencix.com](mailto:sales@sciencix.com) to report the issue. We will assist you with a resolution in a timely fashion.