

# CTS-22008

Deuterium Lamp (2000 hr), i-Series, PDA  
Comparable to: Shimadzu® 228-55626-01



## Technical Data Sheet

---

### 1.0 Description

This Lamp is designed to be compatible with LC-2030, LC-2040, LC-2050 and LC-2060 HPLC detectors. They are comparable to Shimadzu® part number 228-55626-01.

### 2.0 Specifications

#### 2.1 Operating Specifications

Lifetime: 2000 hours

#### 2.2 Materials

Socket Contacts: Gold Plated

Lamp Envelope: UV Glass

### 3.0 Replacement Frequency

Sciencix recommends replacing the lamp when it fails to ignite repeatedly, when the detector fails to calibrate, or when noise and drift exceed application limits.

### 4.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.



### 4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at [sales@sciencix.com](mailto:sales@sciencix.com) to report the issue. We will assist you with a resolution in a timely fashion.